

DATA VALIDATION REPORT ADDENDUM

SEATTLE IRON AND METALS DUST MONITORING - PHASE II

SAMPLING EVENTS JAN-FEB-MARCH 2021 AND APR-MAY-JUNE 2021

Prepared for:

Farallon Consulting 975 5th Ave NW Issaquah, WA 98027

Prepared by:

EcoChem, Inc. 500 Union Street, Suite 1010 Seattle, WA 98101

EcoChem Project: C28201-2

Revised March 29, 2022

Approved for Release:

Christine Ransom Senior Project Chemist EcoChem, Inc.

PROJECT NARRATIVE

Basis for the Data Validation

This report summarizes the results of full validation (EPA Stage 4) performed on PUF cartridge and associated quality control sample data for the Seattle Iron and Metals Dust monitoring project. The report has been revised to remove the "J-20" qualifier from the PCB-004 result for Sample SITE 5 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN). The laboratory "E" flag indicating a result that exceeded the highest calibration standard was in error. A complete list of samples is provided in the **Sample Index**.

Analyses were performed by ALS Life Sciences, Burlington, Ontario. The analytical methods and EcoChem project chemists are listed in the following table:

Analysis	Method	Primary Review	Secondary Review	
Dioxins/Furans	EPA TO-09A	[Clayton	A Rodkin/C Dancom	
PCB Congeners	EPA 1668C	E. Clayton	A. BOOKIN/C. Ransom	

The data were reviewed using guidance and quality control criteria documented in the analytical methods; *Seattle Iron & Metals Dust Monitoring Plan: Phase II* (June 2020); National *Functional Guidelines for Organic Superfund Methods Data Review* (USEPA 2017); and *National Functional Guidelines for High Resolution Superfund Methods Data Review* (USEPA 2016).

EcoChem's goal in assigning data assessment qualifiers is to assist in proper data interpretation. If values are estimated (J or UJ), data may be used for site evaluation and risk assessment purposes but reasons for data qualification should be taken into consideration when interpreting sample concentrations. Data that have been rejected are flagged with (R). Rejected data should not be used for any purpose. If values have no data qualifier assigned, then the data meet the data quality objectives as stated in the documents and methods referenced above.

Validation criteria are included as **Appendix A**. The qualified data summary table (QDST) is included as **Appendix B**. Data Validation Worksheets and project associated communications will be kept on file at EcoChem, Inc. A qualified laboratory electronic data deliverable (EDD) is also submitted.

Sample Index Seattle Iron and Metals - Dust Monitoring Phase II

				РСВ
SDG	Sample ID	Laboratory ID	Dioxins	Congeners
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	\checkmark	\checkmark
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	\checkmark	\checkmark
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	\checkmark	\checkmark
L2569163	SITE 4 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-4	\checkmark	\checkmark
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	\checkmark	\checkmark
L2602390	SITE 1 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-1	\checkmark	\checkmark
L2602390	SITE 2 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-2	\checkmark	\checkmark
L2602390	SITE 3 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-3	\checkmark	\checkmark
L2602390	SITE 4 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-4	\checkmark	\checkmark
L2602390	SITE 5 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-5	\checkmark	\checkmark

DATA VALIDATION REPORT Seattle Iron and Metals - Dust Monitoring Phase II Dioxin/Furan Compounds by EPA TO-09A

This report documents the review of analytical data from the analyses of PUF cartridge samples and the associated laboratory quality control (QC) samples. Samples were analyzed by ALS Life Sciences, Burlington, Ontario. Refer to the **Sample Index** for a complete list of samples.

SDG	Number of Samples and Matrix	VALIDATION LEVEL
L2569163	5 PUF Composites	Stage 4
L2602390	5 PUF Composites	Stage 4

DATA PACKAGE COMPLETENESS

The laboratory submitted all required deliverables. The laboratory followed adequate corrective action processes and all anomalies were discussed in the case narrative.

EDD TO HARDCOPY VERIFICATION

Ten percent (10%) of the results in the laboratory EDD were verified by comparison to the laboratory data package. With the noted exceptions, no errors were noted.

The collection dates and times were missing from the EDD for all samples. Dates were added during validation. The values for the final week of collection were used for the composites.

SDG L2602390: For Sample Site 1-Composite 6 (Wet Season Apr, May, Jun), two corrections were made to the EDD:

- 1,2,3,4,7,9-HxCDF: the value in the result_value column was changed from 0 to 1.3.
- Total HpCDF: the laboratory qualifier was missing from the lab_qualifiers column. A "U" was added.

TECHNICAL DATA VALIDATION

The quality control (QC) requirements that were reviewed are listed in the following table.

1	Sample Receipt, Preservation, and Holding Times	1	Laboratory Control Samples
\checkmark	System Performance and Resolution Checks	1	Matrix Spike/Matrix Spike Duplicates (MS/MSD)
\checkmark	Initial Calibration (ICAL)	1	Field Duplicates
1	Calibration Verification (CCAL)	\checkmark	Reporting Limits and Sample Quantitation
2	Laboratory Blanks	\checkmark	Target Analyte List
1	Field Standard Labeled Compounds	2	Compound Identification
1	Extraction Standard Labeled Compounds	1	Calculation Verification
\checkmark	Cleanup Standards		

✓ Method quality objectives (MQO) and QC criteria have been met. No outliers are noted or discussed. 1 Quality control results are discussed below, but no data were qualified.

2 Quality control outliers that impact the reported data were noted. Data qualifiers were issued as discussed below.

Sample Receipt, Preservation, and Holding Times

The validation guidance documents state that the cooler temperatures should be within an advisory temperature range of $\leq 6^{\circ}$ C.

The validation guidance documents state that the cooler temperatures should be within an advisory temperature range of ≤ 6 °C. Weekly samples were stored frozen until samples for all 4 individual sampling events for the month were collected. These were then composited into a single sample. A holding time of one year was used to evaluate the sample extraction and analysis.

SDG L2569163: Several sample cooler temperatures were greater than the upper control limit, ranging from 2.0-19.7°C. These outliers did not impact data quality; no data were qualified.

SDG L2602390: Several sample cooler temperatures were greater than the upper control limit, ranging from 6.0-15.8°C. These outliers did not impact data quality; no data were qualified.

Calibration Verification

All calibration verification (CCAL) criteria were met. The CCAL were analyzed at the beginning of each sequence as required. All compound concentrations fell within the acceptance limits specified in the method. All ion ratios were acceptable. The S/N ratio was greater than 10, as required, for all compounds. All relative retention times for all target compounds met the required criteria.

Laboratory Blanks

To assess the impact of any blank contaminant on the reported sample results, an action level is established at five times (5x) the concentration reported in the blank. If a contaminant is reported in an associated field sample and the concentration is less than the action level, the result is qualified as not detected (U-7). No action is taken if the sample result is greater than the action level, or for non-detected results. The laboratory assigned R-flags to values when a peak was detected but did not meet identification criteria. These values cannot be considered as positive identifications but are

"estimated maximum possible concentrations" or EMPC. When these occurred in the method blank the results were considered as false positives. No action levels were established for these analytes.

Method blanks (media blanks and reagent blanks) were analyzed at the appropriate frequency. For detected results, the highest concentrations from the two blanks were used to evaluate field sample data. Several target analytes were detected in each blank for every batch. The associated detected results were greater than the 5x action levels, with the following exceptions.

SDG L2602390: Results for 1,2,3,4,7,8-HxCDF in two samples were less than the action level and were qualified as not detected (U-7).

Field Standard Labeled Compounds

Five labeled compounds were added to the sample cartridges prior to sampling. The percent recovery (%R) values method specified control limits are 70-130%. These surrogates are used to evaluate the sampling system and are not used for sample quantitation. No action was taken unless three or more of the compounds were outside of the control limits. No qualifiers were assigned.

Extraction Standard Labeled Compounds

Isotope-stable labeled compounds were added to each field and QC sample. All recoveries were within the method criteria of 40-130% (25-130% for hepta and octa-chlorinated compounds).

Laboratory Control Samples

Laboratory control samples were analyzed at the proper frequency. All recovery values were within the control limits.

Matrix Spikes/Matrix Spike Duplicates

Matrix spike/matrix spike duplicates (MS/MSD) were not analyzed. These are not required by the method. Accuracy was assessed using labeled compound and LCS recoveries. Precision within an analytical batch could not be evaluated.

Field Duplicates

No field duplicates were collected.

Compound Identification

The laboratory reported EMPC or "estimated maximum possible concentrations" values for one or more of the target analytes in all samples. An EMPC value was reported when a peak was detected but did not meet identification criteria as required by the method; therefore, the result cannot be considered as positive identification for the analyte. The lab flagged these results "R". All EMPC results were qualified as not-detected (U-25) at the reported concentrations.

The laboratory uses DB5 MS column, which provides adequate resolution of the TCDF isomers as indicated by the acceptable peak to valley ratios. There were no positive results for 2,3,7,8-TCDF in the field samples. No second column confirmation was necessary.

Calculation Verification

Several results were verified by recalculation from the raw data. No calculation or transcription errors were found.

OVERALL ASSESSMENT

As determined by this evaluation, the laboratory performed the specified analytical method. Accuracy was acceptable as demonstrated by the labeled compound and laboratory control sample recoveries. Precision within an analytical batch could not be evaluated.

Data were qualified as not detected due to method blank contamination. Detection limits were elevated based on ion ratio outliers.

All data, as qualified, are acceptable for use.

DATA VALIDATION REPORT Seattle Iron and Metals - Dust Monitoring Phase II PCB Congeners by EPA 1668C

This report documents the review of analytical data from the analyses of PUF cartridge samples and the associated laboratory quality control (QC) samples. Samples were analyzed by ALS Life Sciences, Burlington, Ontario. Refer to the Sample Index for a complete list of samples.

SDG	Number of Samples and Matrix	Validation Level
L2569163	5 PUF Composites	Stage 4
L2602390	5 PUF Composites	Stage 4

DATA PACKAGE COMPLETENESS

The laboratory submitted all required deliverables. The laboratory followed adequate corrective action processes and all anomalies were discussed in the case narrative.

EDD TO HARDCOPY VERIFICATION

Ten percent (10%) of the results in the laboratory EDD were verified by comparison to the laboratory data package. No errors were noted.

The collection dates and times were missing from the EDD for all samples. Dates were added during validation. The values for the final week of collection were used for the composites.

TECHNICAL DATA VALIDATION

This report documents the review of analytical QC requirements as listed in the following table.

1	Sample Receipt, Preservation, and Holding Times	1	Laboratory Control Samples
\checkmark	Instrument Performance Check	1	Matrix Spike/Matrix Spike Duplicates (MS/MSD)
\checkmark	Initial Calibration (ICAL)	1	Field Duplicates
\checkmark	Continuing Calibration Verification (CCV)	\checkmark	Reporting Limits and Sample Quantitation
2	Laboratory Blanks	\checkmark	Target Analyte List
1	Field Standard Labeled Compounds	2	Compound Identification
1	Extraction Standard Labeled Compounds	\checkmark	Compound Quantitation
\checkmark	Cleanup Standards	1	Calculation Verification

 \checkmark Method quality objectives (MQO) and QC criteria have been met. No outliers are noted or discussed.

1 Quality control results are discussed below, but no data were qualified.

2 Quality control outliers that impact the reported data were noted. Data qualifiers were issued as discussed below.

Sample Receipt, Preservation, and Holding Times

The validation guidance documents state that the cooler temperatures should be within an advisory temperature range of $\leq 6^{\circ}$ C. Weekly samples were stored frozen until 4 weeks-worth of samples

were collected. These were then composited into a single sample. A holding time of one year was used to evaluate the sample extraction and analysis.

SDG L2569163: One sample cooler temperature was greater than the upper control limit at 19.7°C. This outlier did not impact data quality; no data were qualified.

SDG L2602390: Sample cooler temperatures were greater than the upper control limit, ranging from 6.0-15.8°C. These outliers did not impact data quality; no data were qualified.

Laboratory Blanks

To assess the impact of any blank contaminant on the reported sample results, an action level is established at five times (5x) the concentration reported in the blank. If a contaminant is reported in an associated field sample and the concentration is less than the action level, the result is qualified as not detected (U-7). No action is taken if the sample result is greater than the action level, or for non-detected results. The laboratory assigned R-flags to values when a peak was detected but did not meet identification criteria. These values cannot be considered as positive identifications but are "estimated maximum possible concentrations" or EMPC. When these occurred in the method blank, the results were considered as false positives. No action levels were established for these analytes.

Method blanks (media blanks and reagent blanks) were analyzed at the appropriate frequency. For detected results, the highest concentrations from the two blanks were used to evaluate field sample data. Several target analytes were detected in each blank for every batch. All associated detected results were greater than the 5x action levels with the following exceptions.

SDG L2602390: Results for PCB-184 in all samples were less than the action level and were qualified as not detected (U-7).

Field Standards Labeled Compounds

Three labeled compounds were added to the sample cartridges prior to sampling. The percent recovery (%R) control limits are 70-130%. These surrogates are used to evaluate the sampling system. No action was taken if only one of the recoveries was outside of the control limits. If the recoveries indicated a potential problem with the sampling set-up, all results for the sample were estimated (J/UJ-13).

SDG L2569163: The %R values for 13C12-PCB-095 were less than the lower control limit for Samples Site 1-Composite 5 (Wet Season-Jan, Feb, Mar) and Site 3-Composite 5 (Wet Season-Jan, Feb, Mar). The %R values for the other two field standard labeled compounds were acceptable. No qualifiers were assigned.

Extraction Standard Labeled Compounds

Labeled compounds were added as internal standards to each sample prior to extraction. With the exceptions noted below, labeled compound recoveries were within the method specified control limits.

SDG L2569163: The %R values for labeled compounds 13C12-PCB-019 and 13C12-PCB-054 were less than the lower control limit for the method blanks and the laboratory control sample (LCS). Qualifiers are not assigned to QC samples; no action was taken.

Laboratory Control Samples

Laboratory control samples (LCS) were analyzed with each batch. All recoveries were within the control limits of 70%-130%.

Matrix Spikes/Matrix Spike Duplicates

Matrix spike/matrix spike duplicates (MS/MSD) were not analyzed. These are not required by the method. Accuracy was assessed using labeled compound and laboratory control sample recoveries. Precision could not be assessed.

Field Duplicates

No field duplicates were submitted.

Compound Identification

The laboratory reported EMPC or "estimated maximum possible concentrations" values for one or more of the target analytes in all samples. An EMPC value was reported when a peak was detected but did not meet identification criteria as required by the method; therefore, the result cannot be considered as positive identification for the analyte. The lab flagged these results "R". All EMPC results were qualified as not-detected (U-25) at the reported concentrations.

Calculation Verification

Several results were verified by recalculation from the raw data. No transcription or calculation errors were found.

OVERALL ASSESSMENT

As was determined by this evaluation, the laboratory performed the specified analytical method. With the noted exceptions, accuracy was acceptable as demonstrated by the labeled compound and LCS recoveries. Precision within the batches could not be assessed.

Results were qualified as not detected due to method blank contamination. Detection limits were elevated due to ion ratio outliers.

All data, as qualified, are acceptable for use.



APPENDIX A

DATA QUALIFIER DEFINITIONS REASON CODES AND CRITERIA TABLES

Z:\Shared\EC Final Docs\EarthCon 268\26804 Columbus\26804-2\APPENDIX A.docx

DATA VALIDATION QUALIFIER CODES Based on National Functional Guidelines

The following definitions provide brief explanations of the qualifiers assigned to results in the data review process.

U	The analyte was analyzed for, but was not detected above the reported sample quantitation limit.
J	The analyte was positively identified; the associated numerical value is the approximate concentration of the analyte in the sample.
NJ	The analysis indicates the presence of an analyte that has been "tentatively identified" and the associated numerical value represents the approximate concentration.
UJ	The analyte was not detected above the reported sample quantitation limit. However, the reported quantitation limit is approximate and may or may not represent the actual limit of quantitation necessary to accurately and precisely measure the analyte in the sample.
R	The sample results are rejected due to serious deficiencies in the ability to analyze the sample and meet quality control criteria. The presence or absence of the analyte cannot be verified.
The following is an EcoChem	qualifier that may also be assigned during the data review process:

DNR Do not report; a more appropriate result is reported from another analysis or dilution.

DATA QUALIFIER REASON CODES

Group	Code	Reason for Qualification				
Sample Handling	1	Improper Sample Handling or Sample Preservation (i.e., headspace, cooler temperature, pH, summa canister pressure); Exceeded Holding Times				
	24	Instrument Performance (i.e., tune, resolution, retention time window, endrin breakdown, lock-mass)				
	5A	Initial Calibration (RF, %RSD, r ²)				
Instrument Performance	5B	Calibration Verification (CCV, CCAL; RF, %D, %R) Use bias flags (H,L) ¹ where appropriate				
	5C	Initial Calibration Verification (ICV %D, %R) Use bias flags (H,L) ¹ where appropriate				
	6	Field Blank Contamination (Equipment Rinsate, Trip Blank, etc.)				
Blank Contamination	7	Lab Blank Contamination (i.e., method blank, instrument blank, etc.) Use low bias flag (L) ¹ for negative instrument blanks				
	8	Matrix Spike (MS and/or MSD) Recoveries Use bias flags (H,L) ¹ where appropriate				
	9	Precision (all replicates: LCS/LCSD, MS/MSD, Lab Replicate, Field Replicate)				
Precision and Accuracy	10	Laboratory Control Sample Recoveries (a.k.a. Blank Spikes) Use bias flags (H,L) ¹ where appropriate				
	12	Reference Material Use bias flags (H,L) ¹ where appropriate				
	13	Surrogate Spike Recoveries (a.k.a. labeled compounds, recovery standards) Use bias flags (H,L) ¹ where appropriate				
	16	ICP/ICP-MS Serial Dilution Percent Difference				
	17	ICP/ICP-MS Interference Check Standard Recovery Use bias flags (H,L) ¹ where appropriate				
Interferences	19	Internal Standard Performance (i.e., area, retention time, recovery)				
	22	Elevated Detection Limit due to Interference (i.e., chemical and/or matrix)				
	23	Bias from Matrix Interference (i.e. diphenyl ether, PCB/pesticides)				
	2	Chromatographic pattern in sample does not match pattern of calibration standard				
1.1	3	2 nd column confirmation (RPD or %D)				
Quantitation	4	Tentatively Identified Compound (TIC) (associated with NJ only)				
	20	Calibration Range or Linear Range Exceeded				
	25	Compound Identification (i.e., ion ratio, retention time, relative abundance, etc.)				
	11	A more appropriate result is reported (multiple reported analyses i.e., dilutions, re- extractions, etc. Associated with "R" and "DNR" only)				
Miscellaneous	14	Other (See DV report for details)				
	26	Method QC information not provided				

¹H = high bias indicated

L = low bias indicated

QC Element	Acceptance Criteria	Source of Criteria	Action for Non-Conformance		Discussion and Comments
Sample Handling					·
Cooler/Storage Temperature Preservation	Waters/Solids \leq 6°C & in the dark Tissues <-10°C & in the dark Preservation Aqueous: If Cl ₂ is present Thiosulfate must be added and if pH > 9 it must be adjusted to 7 - 9	NFG ⁽¹⁾ Method ⁽²⁾	J(pos)/R(ND) if thiosulfate not added if Cl ₂ present; J(pos)/UJ(ND) if pH not adjusted J(pos)/UJ(ND) if temp > 20°C	1	EcoChem PJ, see TM-05
Holding Time	If properly stored, 1 year or: Extraction (all matrices): 30 days from collection Analysis (all matrices): 45 days from extraction	NFG ⁽¹⁾ Method ⁽²⁾	If not properly stored or HT exceedance: J(pos)/UJ(ND)	1	EcoChem PJ, see TM-05 Gross exceedance = > 1 year 2011 NFG Note: Under CWA, SDWA, and RCRA the HT for H2O is 7 days.
Instrument Performa	nce				
Mass Resolution (Tuning)	PFK (Perfluorokerosene) ≥10,000 resolving power at m/z 304.9824. Exact mass of m/z 380.9760 w/in 5 ppm of theoretical value (380.97410 to 380.97790) . Analyzed prior to ICAL and at the start and end of each 12 hr. shift.	NFG ⁽¹⁾ Method ⁽²⁾	R(pos/ND) all analytes in all samples associated with the tune	24	Notify PM
Windows Defining Mix	Peaks for first and last eluters must be within established retention time windows for each selector group (chlorination level)	NFG ⁽¹⁾ Method ⁽²⁾	If peaks are not completely within windows (clipped): If natives are ok, J(pos)/UJ(ND) homologs (Totals) If natives are affected, R all results for that selector group	24	Notify PM
Column Performance Mix	Both mixes must be analyzed before ICAL and CCAL Valley < 25% (valley = (x/y)*100%) where x = ht. of TCDD (or TCDF) & y = baseline to bottom of valley For all isomers eluting near the 2378-TCDD (TCDF) peak (TCDD only for 8290)	NFG ⁽¹⁾ Method ⁽²⁾	J(pos) if valley > 25%	24	EcoChem PJ, see TM-05, Rev. 2 ; Note: TCDF is evaluated only if second column confirmation is performed
Initial Calibration Sensitivity	S/N ratio > 10 for all native and labeled compounds in CS1 std.	NFG ⁽¹⁾ Method ⁽²⁾	If <10, elevate Det. Limit or R(ND)	5A	
Initial Calibration Selectivity	Ion Abundance ratios within QC limits (Table 8 of method 8290) (Table 9 of method 1613B)	NFG ⁽¹⁾ Method ⁽²⁾	If 2 or more ion ratios are out for one compound in ICAL, J(pos)	5A	EcoChem PJ, see TM-05, Rev. 2

QC Element	Acceptance Criteria	Source of Criteria	Action for Non-Conformance	Reason Code	Discussion and Comments
Instrument Performa	nce (continued)				
Initial Calibration (Minimum 5 stds.) Stability	%RSD < 20% for native compounds %RSD <30% for labeled compounds (%RSD < 35% for labeled compounds under 1613b)	NFG ⁽¹⁾ Method ⁽²⁾	J(pos) natives if %RSD > 20%	5A	
	Absolute RT of ¹³ C ₁₂ -1234-TCDD >25 min on DB5 & >15 min on DB-225	NFG ⁽¹⁾ Method ⁽²⁾	Narrate, no action		EcoChem PJ, see TM-05, Rev. 2
Continuing Calibration (Prior to each 12 hr. shift) Sensitivity	S/N ratio for CS3 standard > 10	NFG ⁽¹⁾ Method ⁽²⁾	If <10, elevate Det. Limit or R(ND)	5B	
Continuing Calibration (Prior to each 12 hr. shift) Selectivity	Ion Abundance ratios within QC limits (Table 8 of method 8290) (Table 9 of method 1613B)	NFG ⁽¹⁾ Method ⁽²⁾	For congener with ion ratio outlier, J(pos) natives in all samples associated with CCAL. No action for labeled congener ion ratio outliers.	25	EcoChem PJ, see TM-05
Continuing Calibration (Prior to each 12 hr. shift)	%D+/-20% for native compounds %D +/-30% for labeled compounds (Must meet limits in Table 6, Method 1613B) If %D in the closing CCAL are within 25%/35%, the mean RF from the two CCAL may be used to calculate samples (Section 8.3.2.4 of 8290).	NFG ⁽¹⁾ Method ⁽²⁾	Labeled compounds: Narrate, no action. Native compounds: 1613: J(pos)/UJ(ND)if %D is outside Table 6 limits J(pos)/R(ND) if %D is +/-75% of Table 6 limits 8290: J(pos)/UJ(ND) if %D = 20% - 75% J(pos)/R(ND) if %D > 75%	5B (H,L) ³	
Stability	Absolute RT of ¹³ C ₁₂ -1234-TCDD and ¹³ C ₁₂ -123789-HxCDD should be ± 15 seconds of ICAL RRT for all other compounds must meet criteria listed in Table 2 Method 1316.	NFG ⁽¹⁾ Method ⁽²⁾	Narrate, no action	5B	EcoChem PJ, see TM-05
Blank Contamination					
Method Blank (MB)	MB: One per matrix per batch of (of ≤ 20 samples) No detected compounds > RL	NFG ⁽¹⁾ Method ⁽²⁾	U(pos) if result is < 5X action level.	7	Hierarchy of blank review: #1 - Review MB, qualify as needed
Field Blank (FB)	FB: frequency as per QAPP No detected compounds > RL	Wethou	U(pos) if result is < 5X action level.	6	#2 - Review FB , qualify as needed

QC Element	Acceptance Criteria	Source of Criteria	Action for Non-Conformance	Reason Code	Discussion and Comments
Precision and Accura	cy				
MS/MSD (recovery)	MS/MSD not typically required for HRMS analyses. If lab analyzes MS/MSD then one set per matrix per batch (of ≤ 20 samples) Use most current laboratory control limits	EcoChem standard policy	J(pos) if both %R > UCL - high bias J(pos)/UJ(ND) if both %R < LCL - low bias J(pos)/R(ND) if both %R < 10% - very low bias J(pos)/UJ(ND) if one > UCL & one < LCL, with no bias PJ if only one %R outlier	8 (H,L) ³	No action if only one spike %R is outside criteria. No action if parent concentration is >4x the amount spiked. Qualify parent sample only unless other QC indicates systematic problems.
MS/MSD (RPD)	MS/MSD not typically required for HRMS analyses. If lab analyzes MS/MSD then one set per matrix per batch (of ≤ 20 samples) Use most current laboratory control limits	EcoChem standard policy	J(pos) in parent sample if RPD > CL	9	Qualify parent sample only.
LCS (or OPR)	One per lab batch (of ≤ 20 samples) Use most current laboratory control limits or Limits from Table 6 of 1613B	NFG ⁽¹⁾ Method ⁽²⁾	J(pos) if %R > UCL - high bias J(pos)/UJ(ND) if %R < LCL - low bias J(pos)/R(ND) if %R < 10% - very low bias	10 (H,L) ³	No action if only one spike %R is outside criteria, when LCSD is analyzed. Qualify all associated samples.
LCS/LCSD (RPD)	LCSD not typically required for HRMS analyses. One set per matrix and batch of 20 samples RPD < 35%	Method ⁽²⁾ Ecochem standard policy	J(pos) assoc. compound in all samples if RPD > CL	9	Qualify all associated samples.
Lab Duplicate (RPD)	Lab Dup not typically required for HRMS analyses. One per lab batch (of ≤ 20 samples) Use most current laboratory control limits	EcoChem standard policy	J(pos)/UJ(ND) if RPD > CL	9	
Labeled Compounds (Internal Standards)	Added to all samples %R = 40% - 135% in all samples 8290 %R must meet limits in Table 7 Method 1613B	NFG ⁽¹⁾ Method ⁽²⁾	J(pos) if %R > UCL - high bias J(pos)/UJ(ND) if %R < LCL - low bias J(pos)/R(ND) if %R < 10% - very low bias	13 (H,L) ³	
Field Duplicates	Solids: RPD <50% OR difference < 2X RL (for results < 5X RL) Aqueous: RPD <35% OR difference < 1X RL (for results < 5X RL)	EcoChem standard policy	Narrate and qualify if required by project	9	Use professional judgment

QC Element	Acceptance Criteria	Source of Criteria	Action for Non-Conformance		Discussion and Comments
Compound ID and Ca	alculation				•
Quantitation/ Identification	All ions for each isomer must maximize within ± 2 seconds. S/N ratio >2.5 Ion ratios must meet criteria listed in Table 8 Method 8290, or Table 9 of 1613B; RRTs w/in limits in Table 2 of 1613B	NFG ⁽¹⁾ Method ⁽²⁾	Narrate in report; qualify if necessary NJ(pos) for retention time outliers. U(pos) for ion ratio outliers.	25	EcoChem PJ, see TM-05
EMPC (estimated maximum possible concentration)	If quantitation identification criteria are not met, laboratory should report an EMPC value.	NFG ⁽¹⁾ Method ⁽²⁾	If laboratory correctly reported an EMPC value, qualify the native compound U(pos) to indicate that the value is a detection limit and qualify total homolog groups J (pos)	25	Use professional judgment See TM-18
Interformers	Interferences from chlorodiphenyl ether compounds	NFG ⁽¹⁾ Method ⁽²⁾	J(pos)/UJ(ND) if present	23	See TM-16
Interferences	Lock masses must not deviate ± 20% from values in Table 8 of 1613B	Method ⁽²⁾	J(pos)/UJ(ND) if present	24	See TM-17
Second Column Confirmation	All 2,3,7,8-TCDF hits must be confirmed on a DB-225 (or equiv) column. All QC criteria must also be met for the confirmation analysis.	NFG ⁽¹⁾ Method ⁽²⁾	Report the DB-225 value. If not performed use PJ.	3	DNR-11 DB5 result if both results from both columns are reported. EcoChem PJ, see TM-05
Calculation Check	Check 10% of field & QC sample results	EcoChem standard policy	Contact laboratory for resolution and/or corrective action	na	Full data validation only.
Electronic Data Deliv	rerable (EDD)			-	
Verification of EDD to hardcopy data	EcoChem verify @ 10% unless problems noted; then increase level up to 100% for next several packages.		Depending on scope of problem, correct at EcoChem (minor issues) to resubmittal by laboratory (major issues).	na	EcoChem Project Manager and/or Database Administrator will work with lab to provide long-term corrective action.
Dilutions, Re- extractions and/or Reanalyses	Report only one result per analyte	Standard reporting policy	Use "DNR" to flag results that will not be reported.	11	

(pos) - positive (detected) results; (ND) - not detected results

¹ National Functional Guidelines for Chlorinated Dibenzo-p-Dioxins (CDDs) & Chlorinated Dibenzofurans (CDFs) Data Review, September 2011

² Polychlorinated Dibenzodioxins (PCDDs) and Polychlorinated Dibenzofurans (PCDFs) by High-Resolution Gas Chromatography/High-Resolution Mass Spectrometry (HRGC/HRMS), USEPA SW-846, Method 8290

² EPA Method 1613, Rev.B, Tetra-through Octa-Chlorinated Dioxins and Furans by Isotope Dilution HRGS/HRMS, October 1994

³ NFG 2013 suggests using "+ / -" to indicate bias; EcoChem has chosen "H" = high bias indicated; "L" = low bias indicated.

PCB Congener Analysis by HRMS (Based on EPA DV Guidance¹ and Method EPA 1668C)

QC Element	Acceptance Criteria	Acceptance Criteria Source of Criteria Action for Non-Conformance			
Sample Handling	•				
Cooler/Storage Temperature Preservation	Waters/Solids ≤ 6°C & in the dark Tissues <-10°C & in the dark Preservation Aqueous: If Cl ₂ is present Thiosulfate must be added and if needed adjust pH to 2 - 3 (drinking water requirement)	EPA ⁽¹⁾ Method ⁽²⁾	J(pos)/R(ND) if thiosulfate not added if Cl ₂ present and J(pos)/UJ(ND) if pH not adjusted; J(pos)/UJ(ND) if temp > 20°C	1	Note: EPA DV guidance documents use < 4°C, method uses ≤ 6°C. Info in EcoChem TM-05 also generally applies.
Holding Time	If properly stored, 1 year prior to extraction. If extracts properly stored (< -10°C & in dark), 1 year from extraction to analysis.	EPA ⁽¹⁾ Method ⁽²⁾	If not properly stored or HT exceeded: J(pos)/UJ(ND)	1	May be dictated by QAPP Info in EcoChem TM-05 also generally applies
Instrument Performa	ince				
Mass Resolution (Tuning)	≥10,000 resolving power at m/z 330.9792 <5 ppm deviation from each m/z listed in Table 7 of method. Analyzed prior to ICAL and at the beginning and end of each 12 hr. shift	EPA ⁽¹⁾ Method ⁽²⁾	R all analytes in all samples associated with a failed tune	24	PFK (Perfluorokerosene) tuning compound
Column Resolution	Mix of all 209 PCBs run prior to each ICAL/12 hours RT of PCB209 must be > 55 min PCB156 & 157 must coelute w/in 2 sec PCB34 & 23 and PCB187 & 182 must be resolved where ((x/y)*100%) < 40% x = ht of valley and y = ht of shortest peak RRT of all congeners must fall within the range in Table 2 of the method	EPA ⁽¹⁾ Method ⁽²⁾	If criteria are not met, review sample chromatograms to determine if sample results are negatively impacted. If so, discuss with client for possible reanalyses, or J(pos) all data.	24	Criteria are for SPB-octyl column. If different column used, see Section 6.9.1.2 of method. Appendix A provides info for DB-1 column
Initial Calibration Sensitivity	S/N ratio > 10 for all native and labeled congeners in CS1 std.	EPA ⁽¹⁾ Method ⁽²⁾	If <10, elevate Det. Limit or R(ND)	5A	
Initial Calibration Selectivity	Ion Abundance ratios within QC limits (Table 8 of Method 1668C)	EPA ⁽¹⁾ Method ⁽²⁾	If ion ratios are out for a given congener in 2 or more standards in ICAL, J(pos) results for that congener in all samples	5A	Professional judgement. The info in EcoChem TM-05 also generally applies
Initial Calibration (Minimum 5 stds.) Stability	%RSD < 20% for congeners listed in Table 3 of method RRT of all congeners must meet Table 2 of method	EPA ⁽¹⁾ Method ⁽²⁾	J(pos) natives if %RSD > 20% RRT outliers: narrate, no action	5A	RRT outliers: professional judgement. The info in EcoChem TM-05 also generally applies
Continuing Calibration (Prior to each 12 hr. shift) Sensitivity	S/N ratio for CS3 standard > 10	EPA ⁽¹⁾ Method ⁽²⁾	If <10, elevate Det. Limit to lowest calibration or R(ND)	5B	
Continuing Calibration (Prior to each 12 hr. shift) Selectivity	Ion Abundance ratios within QC limits (Table 8 of Method 1668C)	EPA ⁽¹⁾ Method ⁽²⁾	No action if %D acceptable, review sample ion ratios, U(pos) if ion ratio outside limits	5B	Professional judgement. The info in EcoChem TM-05 also generally applies.

PCB Congener Analysis by HRMS (Based on EPA DV Guidance¹ and Method EPA 1668C)

QC Element	Acceptance Criteria	Source of Criteria	Action for Non-Conformance	Reason Code	Discussion and Comments	
Continuing Calibration (Prior to each 12 hr. shift)	Recoveries must meet VER% limits in Table 6, Method 1668C	EPA ⁽¹⁾ Method ⁽²⁾	Labeled congeners: Narrate, no action. Native congeners: J(pos)/UJ(ND) for low bias J(pos) for high bias	5B (H,L) ³		
shift) Stability	Absolute RT of all Labeled congeners and Window Defining Congeners must be +/- 15 sec of RT in ICAL RRT of all congeners must be within range in Table 2 of method	EPA ⁽¹⁾ Method ⁽²⁾	Narrate, no action	5B	Professional judgement. The info in EcoChem TM-05 also generally applies	
Blank Contamination						
Method Blank (MB)	MB: One per matrix per batch of (of ≤ 20 samples) No detected congeners	EPA ⁽¹⁾	U(pos) if sample result is < 5X blank concentration	7	Heirarchy of blank review: #1 - Review MB, quaify as needed #2 - Review FB , qualify as needed	
Field Blank (FB)	FB: frequency as per QAPP No detected congeners	Method ⁽²⁾	U(pos) if sample result is < 5X blank concentration	6	EMPC values in blanks as considered to non-detects	
Precision and Accura	cy					
MS/MSD (recovery)	MS/MSD not typically required for HRMS analyses. If lab analyzes MS/MSD then one set per matrix per batch (of ≤ 20 samples) Use most current laboratory control limits	EcoChem standard policy	J(pos) if both %R > UCL - high bias J(pos)/UJ(ND) if both %R < LCL - low bias J(pos)/R(ND) if both %R < 10% - very low bias J(pos)/UJ(ND) if one > UCL & one < LCL, with no bias PJ if only one %R outlier	8 (H,L) ³	No action if only one spike %R is outside criteria. No action if parent concentration is >4x the amount spiked. Qualify parent sample only unless other QC indicates systematic problems.	
MS/MSD (RPD)	MS/MSD not typically required for HRMS analyses. If lab analyzes MS/MSD then one set per matrix per batch (of ≤ 20 samples) Use most current laboratory control limits	EcoChem standard policy	J(pos) in parent sample if RPD > CL	9	Qualify parent sample only.	
LCS (or OPR)	One per lab batch (of ≤ 20 samples) %R must meet limits in Table 6 Method 1668C	EPA ⁽¹⁾ Method ⁽²⁾	J(pos) if %R > UCL - high bias J(pos)/UJ(ND) if %R < LCL - low bias J(pos)/R(ND) if %R < 10% - very low bias	10 (H,L) ³	No action if only one spike %R is outside criteria, when LCSD is analyzed. Qualify all associated samples.	
LCS/LCSD (RPD)	LCS/LCSD not typically required for HRMS analyses. If lab analyzes LCS/LCSD then one set per matrix and batch of 20 samples RPD < 35%	EcoChem standard policy		9	Qualify all associated samples.	
Lab Duplicate (RPD) (if required)	Lab Dup not typically required for HRMS analyses. One per lab batch (of ≤ 20 samples) Use most current laboratory control limits	EcoChem standard policy	J(pos)/UJ(ND) if RPD > CL	9	Optional element. Qualify parent sample only.	

PCB Congener Analysis by HRMS (Based on EPA DV Guidance¹ and Method EPA 1668C)

QC Element	Acceptance Criteria	Source of Criteria	Action for Non-Conformance	Reason Code	Discussion and Comments
Labeled congeners (Internal Standards)	Added to all samples %R must meet limits in Table 6 Method 1668C	EPA ⁽¹⁾ Method ⁽²⁾	J(pos) if %R > UCL - high bias J(pos)/UJ(ND) if %R < LCL - low bias J(pos)/R(ND) if %R <5% - very low bias J(pos)/UJ(ND) if %R between 5-10% for two or more labeled compounds in a substitution group (ie, mono, - di-, trichlorinated)- very low bias	13 (H,L) ³	See next tab for labled congener associations as per Table 2 Method 1668
Field Duplicates	Solids: RPD <50% OR difference < 2X RL (for results < 5X RL) Aqueous: RPD <35% OR difference < 1X RL (for results < 5X RL)	EcoChem standard policy	Narrate and qualify if required by project (EcoChem PJ)	9	RPD values may be dictated by QAPP 35% and 50% are EcoChem defaults
Compound ID and Ca	lculation			•	
Quantitation/ Identification	All ions for each isomer must maximize within +/- 2 seconds. S/N ratio >2.5 Ion ratios must meet criteria listed in Table 8 of 1668C; RRTs w/in limits in Table 2 of 1668C	EPA ⁽¹⁾ Method ⁽²⁾	Narrate in report; qualify if necessary NJ(pos) for retention time outliers. U(pos) for ion ratio outliers.	25	The info in EcoChem TM-05 also generally applies
EMPC (estimated maximum possible concentration)	If quantitation identification criteria are not met, laboratory should report an EMPC value.	EPA ⁽¹⁾ Method ⁽²⁾	If laboratory correctly reported an EMPC value, qualify the native congener U to indicate that the value is an elevated detection limit and qualify total homolog groups J(+)	25	Use professional judgment. See TM-18
Interferences	Lock masses must not deviate +/- 20% from values in Table 7 of 1668C	Method ⁽²⁾	J(pos)/UJ(ND) if present	24	Use professional judgment. See TM-17
Calibration Range	Results greater than highest calibration standard	EcoChem standard policy	Qualify J (pos)	20	If result from dilution analysis is not reported.
Calculation Check	Check 10% of field & QC sample results	EcoChem standard policy	Contact laboratory for resolution and/or corrective action	na	Full data validation only.
Electronic Data Delive	erable (EDD)				
Verification of EDD to hardcopy data	EcoChem verify @ 10% unless problems noted; then increase level up to 100% for next several packages.		Depending on scope of problem, correct at EcoChem (minor issues) to resubmittal by laboratory (major issues).	na	EcoChem Project Manager and/or Database Administrator will work with lab to provide long-term corrective action.
Dilutions, Re-extractions and/or Reanalyses	Report only one result per analyte	Standard reporting policy	Use "DNR" to flag results that will not be reported.	11	

¹ USEPA Region 2 Data Validation, Standard Operating Procedure for EPA Method 1668A, Revision 1, September 2008

(pos): Positive Result(s) (ND): Non-detects

USEPA Region 3 Interim Guidelines for the Validation of Data Generated Using Method 1668 PCB Congener Data, Revision 0, April 2004 USEPA Region 10 SOP For the Validation of Method 1668 Toxic, Dioxin-like, PCB Data, Revision 1, December 1995

² EPA Method 1668, Rev.C, Chlorinated Biphenyl Congeners in Water, Soil, Sediment, Biosolids, and Tissue by HRGC/HRMS, April 2010

³ "H" = high bias indicated; "L" = low bias indicated

PCB by 1668C Labeled Compound

																156L/									
1	L 3L	4L	15L	19L	37L	54L	77L	81L	104L	105L	114L	118L	123L	126L	155L	157L	167L	169L	188L	189L	202L	205L	206L	208L	209L
	I 2	4	5	16	16	40	40	40	82	82	82	82	82	126	128	128	128	128	170	170	194	194	206	207	209
2	2 3	5	6	17	17	41	41	41	83	83	83	83	83		129	129	129	129	171	171	195	195	207	208	
		6	7	18	18	42	42	42	84	84	84	84	84		130	130	130	130	172	172	196	196			
		7	8	19	20	43	43	43	85	85	85	85	85		131	131	131	131	173	173	197	197			
		8	9	20	21	44	44	44	86	86	86	86	86		132	132	132	132	174	174	198	198			
		9	10	21	22	45	45	45	87	87	87	87	87		133	133	133	133	175	175	199	199			
		10	11	22	23	46	46	46	88	88	88	88	88		134	134	134	134	176	176	200	200			
		11	12	23	24	47	47	47	89	89	89	89	89		135	135	135	135	177	177	201	201			
		12	13	24	25	48	48	48	90	90	90	90	90		136	136	136	136	178	178	202	203			
		13	14	25	26	49	49	49	91	91	91	91	91		137	137	137	137	179	179	203	204			
		14	15	26	27	50	50	50	92	92	92	92	92		138	138	138	138	180	180	204	205			
				27	28	51	51	51	93	93	93	93	93		139	139	139	139	181	181					
				28	29	52	52	52	94	94	94	94	94		140	140	140	140	182	182					
				29	30	53	53	53	95	95	95	95	95		141	141	141	141	183	183					
				30	31	54	55	55	96	96	96	96	96		142	142	142	142	184	184					
				31	32	55	56	56	97	97	97	97	97		143	143	143	143	185	185					
				32	33	56	57	57	98	98	98	98	98		144	144	144	144	186	186					
				33	34	57	58	58	99	99	99	99	99		145	145	145	145	187	187					
				34	35	58	59	59	100	100	100	100	100		146	146	146	146	188	189					
				35	36	59	60	60	101	101	101	101	101		147	147	147	147	190	190					
				30	31	60	61	61	102	102	102	102	102		148	148	148	148	191	191					
				38	38	61	62	62	103	103	103	103	103		149	149	149	149	192	192					
				39	39	62	63	63	104	105	100	100	100		150	150	150	150	193	193]				
						64	65	65	106	100	107	107	107		151	151	151	151							
						65	66	66	107	107	100	100	100		152	152	152	152							
						66	67	67	100	100	1109	109	109		153	153	153	153							
						67	68	68	110	110	111	111	111		155	154	158	158							
						68	60	60	111	111	112	112	112		158	157	150	150	•						
						69	70	70	112	112	113	113	113		159	158	160	160	1						
						70	71	71	113	113	114	115	115		160	159	161	161							
						71	72	72	115	115	115	116	116		161	160	162	162							
						72	73	73	116	116	116	117	117		162	161	163	163	1						
						73	74	74	117	117	117	118	119		163	162	164	164	1						
						74	75	75	119	119	119	119	120		164	163	165	165							
						75	76	76	120	120	120	120	121	1	165	164	166	166	1						
						76	77	78	121	121	121	121	122	1	166	165	167	168	1						
						78	78	79	122	122	122	122	123	1	168	166	168	169	1						
						79	79	80	124	124	124	124	124	1		168		-	•						
						80	80	81	125	125	125	125	125	1		P	-								
						8	-	-	127	127	127	127	127	1											



APPENDIX B

QUALIFIED DATA SUMMARY TABLE

								DV	Reason
SDG	Sample ID	Lab ID	Method	Analyte	Result	Units	Lab Flag	Qualifier	Code
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA1668C	PCB-126	22	pg	M,J,R	U	25
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA1668C	PCB-145	5.8	pg	M,J,R	U	25
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA1668C	PCB-148	6.6	pg	M,J,R	U	25
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA1668C	PCB-169	4.5	pg	M,J,R	U	25
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA1668C	PCB-182	5.3	pg	M,J,R	U	25
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA1668C	PCB-185	27	pg	J,R	U	25
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA1668C	PCB-192	4.3	pg	J,R	U	25
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA1668C	PCB-196	47	pg	J,R	U	25
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA1668C	PCB-104	46	pg	M,J,R	U	25
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA1668C	PCB-126	51	pg	M,J,R	U	25
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA1668C	PCB-142	11	pg	J,R	U	25
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA1668C	PCB-169	14	pg	M,J,R	U	25
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA1668C	PCB-186	2	pg	J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA1668C	PCB-068	110	pg	R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA1668C	PCB-104	9.8	pg	J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA1668C	PCB-126	18	pg	M,J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA1668C	PCB-145	8	pg	M,J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA1668C	PCB-169	4.5	pg	J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA1668C	PCB-188	3.3	pg	J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA1668C	PCB-181	8.2	pg	J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA1668C	PCB-189	19	pg	M,J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA1668C	PCB-197	9.8	pg	M,J,R	U	25
L2569163	SITE 4 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-4	EPA1668C	PCB-126	26	pg	M,J,R	U	25
L2569163	SITE 4 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-4	EPA1668C	PCB-169	6	pg	J,R	U	25
L2569163	SITE 4 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-4	EPA1668C	PCB-188	3.7	pg	J,R	U	25
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA-TO-09A	1,2,3,4,7,8-HxCDD	1.8	pg	M,J,R	U	25
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA-TO-09A	1,2,3,6,7,8-HxCDD	3	pg	J,R	U	25
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA-TO-09A	1,2,3,7,8,9-HxCDD	3.1	pg	M,J,R	U	25
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA-TO-09A	2,3,4,7,8-PeCDF	4.1	pg	M,J,R	U	25
L2569163	SITE 1 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-1	EPA-TO-09A	2,3,4,6,7,8-HxCDF	3.1	pg	M,J,R	U	25
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA-TO-09A	1,2,3,7,8-PeCDD	1.9	pg	M,J,R	U	25

								DV	Reason
SDG	Sample ID	Lab ID	Method	Analyte	Result	Units	Lab Flag	Qualifier	Code
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA-TO-09A	1,2,3,6,7,8-HxCDD	4.3	pg	M,J,R	U	25
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA-TO-09A	1,2,3,7,8,9-HxCDD	3.4	pg	M,J,R	U	25
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA-TO-09A	2,3,7,8-TCDF	3.1	pg	J,R	U	25
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA-TO-09A	1,2,3,7,8-PeCDF	3.4	pg	M,J,R	U	25
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA-TO-09A	1,2,3,6,7,8-HxCDF	4.5	pg	M,J,R	U	25
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA-TO-09A	1,2,3,7,8,9-HxCDF	2.4	pg	M,J,R	U	25
L2569163	SITE 2 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-2	EPA-TO-09A	1,2,3,4,7,8,9-HpCDF	0	pg	M,J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA-TO-09A	1,2,3,6,7,8-HxCDD	3	pg	M,J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA-TO-09A	2,3,4,7,8-PeCDF	3	pg	M,J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA-TO-09A	1,2,3,4,7,8-HxCDF	2.6	pg	J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA-TO-09A	2,3,4,6,7,8-HxCDF	3.6	pg	J,R	U	25
L2569163	SITE 3 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-3	EPA-TO-09A	1,2,3,4,6,7,8-HpCDF	10	pg	M,J,R	U	25
L2569163	SITE 4 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-4	EPA-TO-09A	1,2,3,7,8-PeCDF	3.5	pg	J,R	U	25
L2569163	SITE 4 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-4	EPA-TO-09A	1,2,3,4,6,7,8-HpCDF	11	pg	M,J,R	U	25
L2569163	SITE 4 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-4	EPA-TO-09A	OCDF	25	pg	M,J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA-TO-09A	1,2,3,4,7,8-HxCDD	2.5	pg	M,J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA-TO-09A	1,2,3,7,8,9-HxCDD	4.1	pg	M,J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA-TO-09A	1,2,3,7,8-PeCDF	3	pg	M,J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA-TO-09A	1,2,3,4,7,8-HxCDF	3.2	pg	M,J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA-TO-09A	1,2,3,6,7,8-HxCDF	1.9	pg	J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA-TO-09A	1,2,3,4,6,7,8-HpCDF	9.9	pg	J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA-TO-09A	1,2,3,4,7,8,9-HpCDF	1.8	pg	J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA1668C	PCB-121	3.3	pg	J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA1668C	PCB-126	28	pg	M,J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA1668C	PCB-169	4.9	pg	J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA1668C	PCB-184	5.5	pg	J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA1668C	PCB-182	7.2	pg	M,J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA1668C	PCB-181	11	pg	M,J,R	U	25
L2569163	SITE 5 - COMPOSITE 5 (WET SEASON - JAN, FEB, MAR)	L2569163-5	EPA1668C	PCB-205	7.9	pg	M,J,R	U	25
L2602390	SITE 1 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-1	EPA-TO-09A	2,3,4,7,8-PeCDF	2.5	pg	J,R	U	25
L2602390	SITE 1 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-1	EPA-TO-09A	1,2,3,4,7,8-HxCDF	1.3	pg	J,R	U	25

								DV	Reason
SDG	Sample ID	Lab ID	Method	Analyte	Result	Units	Lab Flag	Qualifier	Code
L2602390	SITE 1 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-1	EPA-TO-09A	1,2,3,4,6,7,8-HpCDF	6.8	pg	J,R	U	25
L2602390	SITE 2 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-2	EPA-TO-09A	1,2,3,7,8-PeCDD	1.7	pg	J,R	U	25
L2602390	SITE 2 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-2	EPA-TO-09A	1,2,3,6,7,8-HxCDD	2.1	pg	J,R	U	25
L2602390	SITE 2 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-2	EPA-TO-09A	1,2,3,7,8,9-HxCDD	2.9	pg	J,R	U	25
L2602390	SITE 2 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-2	EPA-TO-09A	2,3,7,8-TCDF	4	pg	R	U	25
L2602390	SITE 2 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-2	EPA-TO-09A	1,2,3,4,7,8-HxCDF	4.8	pg	M,J,R	U	25
L2602390	SITE 2 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-2	EPA-TO-09A	1,2,3,4,6,7,8-HpCDF	23	pg	M,R	U	25
L2602390	SITE 2 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-2	EPA-TO-09A	OCDF	19	pg	J,R	U	25
L2602390	SITE 3 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-3	EPA-TO-09A	1,2,3,4,7,8-HxCDD	2	pg	M,J,R	U	25
L2602390	SITE 3 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-3	EPA-TO-09A	2,3,4,7,8-PeCDF	2.6	pg	M,J,R	U	25
L2602390	SITE 3 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-3	EPA-TO-09A	1,2,3,4,7,8-HxCDF	3.4	pg	M,J,R	U	25
L2602390	SITE 3 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-3	EPA-TO-09A	1,2,3,6,7,8-HxCDF	2.5	pg	M,J,R	U	25
L2602390	SITE 3 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-3	EPA-TO-09A	2,3,4,6,7,8-HxCDF	2	pg	M,J,R	U	25
L2602390	SITE 4 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-4	EPA-TO-09A	2,3,4,7,8-PeCDF	3	pg	M,J,R	U	25
L2602390	SITE 4 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-4	EPA-TO-09A	1,2,3,4,7,8-HxCDF	3.13	pg	J,B	U	7
L2602390	SITE 4 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-4	EPA-TO-09A	1,2,3,6,7,8-HxCDF	1.9	pg	M,J,R	U	25
L2602390	SITE 4 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-4	EPA-TO-09A	1,2,3,4,6,7,8-HpCDF	14	pg	M,J,R	U	25
L2602390	SITE 5 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-5	EPA-TO-09A	1,2,3,4,7,8-HxCDD	1.4	pg	M,J,R	U	25
L2602390	SITE 5 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-5	EPA-TO-09A	2,3,4,7,8-PeCDF	2.8	pg	J,R	U	25
L2602390	SITE 5 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-5	EPA-TO-09A	1,2,3,4,7,8-HxCDF	3.17	pg	M,J,B	U	7
L2602390	SITE 5 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-5	EPA-TO-09A	1,2,3,6,7,8-HxCDF	2.8	pg	J,R	U	25
L2602390	SITE 5 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-5	EPA-TO-09A	1,2,3,4,7,8,9-HpCDF	1.6	pg	J,R	U	25
L2602390	SITE 1 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-1	EPA1668C	PCB-127	5.3	pg	M,J,R	U	25
L2602390	SITE 1 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-1	EPA1668C	PCB-184	6.19	pg	J	U	7
L2602390	SITE 1 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-1	EPA1668C	PCB-181	5.5	pg	J,R	U	25
L2602390	SITE 2 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-2	EPA1668C	PCB-165	12	pg	J,R	U	25
L2602390	SITE 2 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-2	EPA1668C	PCB-184	8.24	pg	M,J	U	7
L2602390	SITE 2 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-2	EPA1668C	PCB-181	30	pg	J,R	U	25
L2602390	SITE 3 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-3	EPA1668C	PCB-078	16	pg	M,J,R	U	25
L2602390	SITE 3 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-3	EPA1668C	PCB-162	45	pg	M,J,R	U	25
L2602390	SITE 3 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-3	EPA1668C	PCB-184	13.3	pg	J	U	7

								DV	Reason
SDG	Sample ID	Lab ID	Method	Analyte	Result	Units	Lab Flag	Qualifier	Code
L2602390	SITE 4 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-4	EPA1668C	PCB-148	25	pg	M,J,R	U	25
L2602390	SITE 4 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-4	EPA1668C	PCB-165	10	pg	J,R	U	25
L2602390	SITE 4 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-4	EPA1668C	PCB-162	55	pg	J,R	U	25
L2602390	SITE 4 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-4	EPA1668C	PCB-184	10.4	pg	J	U	7
L2602390	SITE 5 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-5	EPA1668C	PCB-184	18.2	pg	M,J	U	7
L2602390	SITE 5 - COMPOSITE 6 (WET SEASON - APR, MAY, JUN)	L2602390-5	EPA1668C	PCB-186	2.7	pg	J,R	U	25